

PCR Product Clean-up Protocol

PCR product clean-up in 96 well format is fully automated with our in house liquid handling workstation, epMOTION 5070.

(<http://www.eppendorf.com/epmotion/start.php?nav=51&lang=0>).

We follow Beckman Coulter's Ampure PCR product clean-up kit and protocol.

[http://www.beckmangenomics.com/documents/products/ampure/AMPure_Protocol_Supplement_\(LVR\).pdf](http://www.beckmangenomics.com/documents/products/ampure/AMPure_Protocol_Supplement_(LVR).pdf)